

# PRAVEZA<sup>®</sup> URETERAL ACCESS SHEATH

Used to establish a conduit during urological procedures facilitating the passage for Ureteroscope / Endoscope and other surgical instruments into the urinary tract (RIRS).

## Key Features

- Kink resistant – reinforced with a stainless steel coil designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement of sheath under fluoroscopy
- Sizes available for all types of flexible scopes (Storz, Wolf & Single use flexible scopes)
- Seamless transition between dilator and sheath for easy and atraumatic insertion

## Sizes available for all types of flexible scopes

Storz, Wolf & Single use flexible scopes

## Size (in Fr)

ID	9.5	10	11.5	12
OD	11.5	12	13.5	14

## Length (in cm)

28, 35, 45 & 55

