

N QUADRA[®] Z - NITINOL STONE BASKET ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

2.5, 3, 4 & 5

Length (in cm)

70, 90, 110 & 120



N QUADRA[®] - NITINOL STONE BASKET WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip

Size (in Fr)

2,5, 3, 4 & 5

Length (in cm)

70, 90, 110 & 120



N QUADRA[®] HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

2.5, 3, 4 & 5

Length (in cm)

70, 90, 110 & 120



N QUADRA[®] HELICAL - NITINOL STONE BASKET HELICAL WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

2.5, 3, 4 & 5

Length (in cm)

70, 90, 110 & 120



N QUADRA® FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

Flexible Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx.

Key Features

- Flexible design shaft which can angulate up to 270 degrees as per the flexible scope direction during RIRS Procedure
- Kink resistant braided shaft provides durability, making removal of several fragments possible
- Tipless design helps in easy capturing of stones by allowing the basket to open at the base of calyx and encircle the stone
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

1.8, 1.9, 2, 2.2 & 2.4

Length (in cm)

120

