


P Flex Ureteral Stent (Double J Stent)



Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- One End Multiloop
- With or Without Hydrophilic coating
- With or Without String

Package Includes

Stent with pusher

Also available in set, which includes

Stent with pusher & guidewire

Intended Use:

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features:

- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter & smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization

Standard	
Size (in Fr)	Length (in cm)
4.5, 4.7, 4.8, 5, 5.5, 6, 6.5, 7 & 8	12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Multiloop	
Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	18-24 & 24-30

Pediatric	
Size (in Fr)	Length (in cm)
2.8, 3, 3.5, 3.7, 3.8 & 4	6, 8, 10, 12, 14 & 16