



karex

Foley Balloon Catheters



Foley Balloon Catheters



Manufactured under controlled and clean environment in accordance to ISO 13485 G.M.P. Norm and CE-Mark, all our foley balloon catheters are made from the highest quality of medical grade latex. They are specially designed for transurethral drainage of the urinary bladder, commonly used on patients who are anesthetized or sedated for surgery or other medical care.

We produce a wide range of catheter products which are bio compatible, designed to meet all user and patient requirements.

WHO IS KAREX?

- ◆ *Karex has the market reputation as proven track record as the world's largest condom maker*
- ◆ *We have the capacity to produce over 5 billion condoms annually*
- ◆ *We export to over 130 countries globally*
- ◆ *We are compliant with various management systems such as ISO 13485, ISO 14001, US FDA 21 CFR Part 820, Brazil RDC 16/2013 and GMP*
- ◆ *We also manufacture other products such as condoms, personal lubricants and probe covers*



COMPREHENSIVE RANGE OF CATHETERS

- ◆ Comprehensive range of intermittent, indwelling short term and long term, irrigation, post operative catheters.
- ◆ Made from biocompatible latex with siliconised finishing or silicone elastomer coating.
- ◆ Drainage eyes and tip configuration designed to meet medical profession requirements.
- ◆ Size differentiation by color code and choice of inflation valves.
- ◆ Inclusive range of catheter types, sizes and balloons.




SPECIFIC FEATURES

- ◆ Smooth, round and firm tip to facilitate insertion.
- ◆ Rounded drainage eyes to improve urine flow and minimize encrustation.
- ◆ Variety of catheter colors in translucent, blue, green and red to serve user preferences.
- ◆ Manufacturing code printed on every catheter produced for traceability.
- ◆ Ideal catheter funnel dimension enables connection to fit all urine drainage systems.
- ◆ Short tip and symmetrical balloon shape to minimize bladder irritation.



POST OPERATIVE REINFORCED CATHETER

- ◆ Catheter with re-inforced wall to prevent the collapse of drainage lumen.
- ◆ Tip configuration includes Nelaton, Couvalaire, Dufour, Mercier and Frohmuller designed to meet medical profession requirements.





CATHETER SYSTEM



- ◆ **One Way Catheter**
Red urethral intermittent catheter with stiff shaft to ease insertion.



- ◆ **Two Way Catheter**
Short term indwelling catheter with siliconised finishing to enhance patient comfort.








- ◆ **Three Way Catheter**
3 way siliconised or silicone elastomer coated catheter for irrigation drainage.











- ◆ **Long Term Catheter**
Catheter shaft and drainage lumen are coated with silicone elastomer.







Siliconised Catheter

Type	Tip	Valve	Code	Size: Charriere/French	Balloon Capacity	Length
1 Way 	Nelaton 2 opposed eyes, concave tip 	Standard	110	08/10/12/14/16/ 18/20/22/24/26	00	40cm
2 Way 	Nelaton 2 opposed eyes, concave tip 	Standard Hard Valve Soft Valve	211 212	12/14/16/18/ 20/22/24/26/28/30	10ml & 30ml	40cm
		Female Hard Valve Soft Valve	411 412	12/14/16/18/ 20/22/24/26	10ml & 30ml	22cm
		Paediatric Hard Valve Soft Valve	511 512	06 08/10	3ml 5ml	30cm
	Tiemann 2 opposed eyes, curved, concave tip 	Hard Valve	611	12/14/16/18/ 20/22/24/26	10ml & 30ml	40cm

Siliconised Catheter







Type	Tip	Valve	Code	Size: Charriere/French	Balloon Capacity	Length
3 Way, Irrigation 	Nelaton 2 opposed eyes, concave tip 	Hard Valve Soft Valve	311 312	16/18/20/22/24/26	10ml & 30ml	40cm
2 Way, Reinforced 	Haematuria Nelaton 4 opposed eyes, concave tip 	Hard Valve	221	18/20/22/24/26	30ml	40cm
	Haematuria Couvelaire 4 opposed eyes, "flute" tip 	Hard Valve	231	18/20/22/24/26	30ml	40cm
	Tamponade Dufour 4 opposed eyes, concave "flute" tip 	Hard Valve	241	18/20/22/24/26	60ml	40cm
	Tamponade Mercier 4 opposed eyes, concave tip 	Hard Valve	251	18/20/22/24/26	60ml	40cm
	Tamponade Frohmuller 2 staggered eyes, concave tip 	Hard Valve	261	18/20/22/24/26	60ml	40cm

Silicone Elastomer Coated Catheter, fitted with hard plastic valve

Type	Tip	Colour	Code	Size: Charriere/French	Balloon Capacity	Length
2 Way 	Nelaton 2 opposed eyes, concave tip 	Standard Green Blue	213 217	12/14/16/18/20/ 22/24/26/28	10ml & 30ml	40cm
		Paediatric Green Blue	513 517	06 08/10	3ml 5ml	30cm
	Suprapubic 2 opposed eyes, 1 eye on tip 	Green Blue	283 287	12/14/16/18/ 20/22/24/26	10ml	40cm
	Tiemann 2 opposed eyes, curved, concave tip 	Green Blue	613 617	12/14/16/18/ 20/22/24/26	10ml & 30ml	40cm
3 Way Irrigation 	Nelaton 2 opposed eyes, concave tip 	Green Blue	313 317	16/18/20/22/24/26	10ml & 30ml	40cm



Siliconised Catheter

Type	Tip	Valve	Code	Size: Charriere/French	Balloon Capacity	Length
3 Way, Reinforced 	Haematuria Nelaton 4 opposed eyes, concave tip 	Hard Valve	321	18/20/22/24/26	30ml	40cm
	Haematuria Couvelaire 4 opposed eyes, "flute" tip 	Hard Valve	331	18/20/22/24/26	30ml	40cm
	Tamponade Dufour 4 opposed eyes, concave "flute" tip 	Hard Valve	341	18/20/22/24/26	60ml	40cm
	Tamponade Mercier 4 opposed eyes, concave tip 	Hard Valve	351	18/20/22/24/26	60ml	40cm
	Tamponade Frohmuller 2 staggered eyes, concave tip 	Hard Valve	361	18/20/22/24/26	60ml	40cm



HOW TO ORDER

Code

2 1 1

Size

1 6

Balloon Capacity

1 0

= Catalogue No

= 2111610

Get in touch with us at enquiry@karex.com.my

Marketing Office :
Lot 594, Persiaran Raja Lumu, Pandamaran Industrial Estate,
42000 Port Klang, Selangor, Malaysia
Tel: +603 3165 6688 Fax: +603 3166 2000

karex