

M/L-10

MULTI-FIRE LAPAROSCOPIC CLIP APPLIER



Get a handle on
cost-conscious
sustainability.



M/L-10 MULTI-FIRE LAPAROSCOPIC CLIP APPLIER

The M/L-10 Clip Applier is a reusable multi-fire device with a high-quality reusable handle. It employs disposable cartridges that contain 19 implant-grade titanium clips, the last two of which are color-coded.



LIGATING CLIP CARTRIDGE

- Cross pattern design
- 2 to 5 mm vessel/cystic duct ligation capability
- Surgical-grade titanium
- Tips-first closure sequence for maximum ligating capability and dependability
- Consistent, uniform closure

ORDERING INFORMATION :

M/L-10 Multi-Fire Clip Applier	REF 1002
M/L-10 Clip Cartridge (19 clips)	REF 1112

Cost Saving Procedure Kit:

Clip Cartridge with 3142 Endo Cut Scissor Tip	REF 1242
---	----------



Microline Surgical Inc.
50 Dunham Road Suite 1500 | Beverly, MA 01915 USA
Phone 978-922-9810 | Fax 978-922-9209

Microline EMEA
Storkstraat 17B | 3833 LB Leusden, The Netherlands
Phone +31 (0) 33 76 000 11 | Fax +31 (0) 33 76 000 14

www.microlinesurgical.com

© 2022 Microline Surgical. Microline is a trademark and/or registered trademark of Microline Surgical.
P/N: 09-39-00054-00-US Rev G, 10/22

Important: Prior to use, refer to Instructions for Use supplied with this device for indications, contraindications, adverse events, suggested procedures, warnings and precautions for full product information.

Caution: Federal (USA) law restricts this device to sale by or on the order of a Physician.